

## L2

# Monitoring cell membrane electroporation with ratiometric fluorescent dye Fura-2AM

**Gorazd Pucihar**

*University of Ljubljana, Faculty of Electrical Engineering*

Duration of the experiments: 60 min

Max. number of participants: 4

Location: Cell Culture Laboratory 1

Level: Advanced

### PREREQUISITES

Participants should be familiar with Laboratory safety (S1). The basic knowledge of handling with cells is required for this laboratory practice.

**The aim** of this laboratory practice is to observe electroporation with fluorescent dye and to determine the effect of cell size, shape and orientation on the efficiency of electroporation.


### THEORETICAL BACKGROUND

Exposure of a biological cell to a sufficiently strong external electric field leads to a detectable increase of membrane permeability, a phenomenon termed electropermeabilization. Because it is assumed that an increased permeability is related to the occurrence of hydrophilic pores in the membrane, the phenomenon is often termed electroporation. Provided that the parameters of the electric field (amplitude, duration, number of pulses, frequency) are moderate, the increased permeability is reversible, and cells recover within a few minutes after the exposure. During the state of high permeability the molecules for which the membrane is otherwise impermeable (e.g. drugs, DNA) can be transported across the membrane. Electroporation is nowadays used in biochemistry, molecular biology, and different fields of medicine. It has already become an established method in oncology for electrochemotherapy of tumours, and holds great promises in gene therapy.

The efficiency of electroporation is influenced by the parameters of the electric field, cell size and geometry, and physiological characteristics of the medium surrounding the cell. Different fluorescent dyes (e.g. Propidium Iodide, Lucifer Yellow, Fura-2, Fura-3,...) can be employed to investigate the influence of these parameters on electroporation and the same dyes can be used to monitor electroporation [1, 2, 3, 4].

### EXPERIMENT


We will monitor cell membrane electroporation using a fluorescent calcium sensitive indicator Fura-2AM. Calcium ions, present in the extracellular medium, do not readily cross an intact (nonporated) cell membrane and the intracellular  $\text{Ca}^{2+}$  concentration is low. Once the membrane becomes permeable due to electroporation,  $\text{Ca}^{2+}$  ions enter the cells, where they bind to the dye and change its excitation and emission spectrum (Figure 1) [5].



**Figure 1:** (A) Fluorescence excitation spectra of Fura-2 for different concentrations of  $\text{Ca}^{2+}$  (image from <http://probes.invitrogen.com/handbook/figures/0554.html>). (B) Schematic of the experiment.

**Protocol:** The experiments will be performed on H9c2 rat cardiac myoblast cell line grown in Lab-Tek chambers (TermoFisher Scientific) in Dulbecco's Modified Eagle Medium DMEM supplemented with 10% foetal bovine serum, L-glutamine (all three from Sigma-Aldrich) and antibiotics. Plate  $8 \times 10^4$  cells on cover glass of Lab-Tek chamber and keep them for 24 hours in the incubator. Before experiments, replace the culture medium with fresh medium containing 2  $\mu\text{M}$  Fura-2 AM (TermoFisher Scientific). After 25-30 minutes of incubation at  $37^\circ\text{C}$  wash the excess dye and leave 1.5 ml of culture medium in the chamber. Place the chamber under a fluorescence microscope (Thunder Imager Live Cell system for fluorescence microscopy, Leica Microsystems) and use  $\times 63$  objective. Insert two parallel Pt/Ir wire electrodes with a 5 mm distance between them to the bottom of the chamber. Acquire a pair of two images (with 340, 380 nm excitation wavelength) each second for one minute in time-lapse acquisition mode, using sCMOS Leica fast camera and LASX software. Apply the electric pulse at the 5<sup>th</sup> s of recording: using a BetaTech device, deliver one electric pulse of 100  $\mu\text{s}$  with voltages varying from 150 to 300 V. In LASX software, observe how the ratio image ( $R = F_{340}/F_{380}$ ) changes in time in each cell [3, 4].

Wait for 5 minutes (resealing and recovery) and apply pulse with a higher amplitude. After each pulse, determine which cells are being electroporated (Figure 2). Observe, which cells become electroporated at lower and which at higher pulse amplitudes.



**Figure 2:** Cells stained with Fura-2AM and exposed to electric pulse with increasing amplitude.

## REFERENCES:

1. Valič B., Golzio M., Pavlin M., Schatz A., Faurie C., Gabriel B., Teissié J., Rols M.P., Miklavčič D. Effect of electric field induced transmembrane potential on spheroidal cells: theory and experiment. *Eur Biophys J*, 32:519-528, 2003.
2. Towhidi L., Kotnik T., Pucihar G., Firoozabadi S.M.P., Mozdarani H., Miklavčič D. Variability of the minimal transmembrane voltage resulting in detectable membrane electroporation. *Electromagn Biol Med*, 27:372-385, 2008.
3. Pucihar G., Krmelj J., Reberšek M., Batista Napotnik T., Miklavčič D. Equivalent pulse parameters for electroporation. *IEEE T Biomed Eng*, 58:3279-3288, 2011.
4. Dermol-Černe J., Batista Napotnik T., Reberšek M., Miklavčič D. Short microsecond pulses achieve homogeneous electroporation of elongated biological cells irrespective of their orientation in electric field. *Sci Rep*, 10:9149: 1-17, 2020.

5. Grynkiewicz G., Poenie M., Tsien R.Y. A new generation of Ca<sup>2+</sup> indicators with greatly improved fluorescence properties. *J Biol Chem*, 260:3440-3450, 1985.

---

## NOTES & RESULTS

---